

**Anti-RIPK1 antibody (620-670) {Biotin} (STJ502799)**

STJ502799

**GENERAL INFORMATION**

|                          |   |
|--------------------------|---|
| <b>Product Type</b>      | Primary antibodies  |
| <b>Short Description</b> | Rabbit polyclonal antibody anti-RIPK1 (620-670) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications. |
| <b>Applications</b>      | ELISA, IP, WB   |
| <b>Host/Source</b>       | Rabbit  |
| <b>Reactivity</b>        | Human   |

**PRODUCT PROPERTIES**

|                            |   |
|----------------------------|---|
| <b>Clonality</b>           | Polyclonal  |
| <b>Clone ID</b>            |   |
| <b>Concentration</b>       | 0.55 µg/µl  |
| <b>Conjugation</b>         | Biotin  |
| <b>Purification</b>        | Affinity Purified   |
| <b>Dilution Range</b>      | WB: 1:500<br>DB: 1:10, 000<br>ELISA: 1:10, 000<br>IP: 1:150     |
| <b>Formulation</b>         |   |
| <b>Isotype</b>             | IgG   |
| <b>Storage Instruction</b> | Store at -20°C for long term storage. Avoid freeze-thaw cycles. |

**TARGET INFORMATION**

|                           |   |
|---------------------------|---|
| <b>Gene ID</b>            | <a href="#">8737</a>  |
| <b>Gene Symbol</b>        | <a href="#">RIPK1</a>   |
| <b>Uniprot ID</b>         | <a href="#">RIPK1_HUMAN</a>   |
| <b>Immunogen</b>          | Synthetic peptide taken within amino acid region 620-670 on human Receptor-interacting Serine/Threonine-protein kinase 1 protein. |
| <b>Immunogen Region</b>   | 620-670   |
| <b>Specificity</b>        |   |
| <b>Immunogen Sequence</b> |   |