

Anti-RIPK2 antibody (410-460) {Biotin} (STJ502802)

STJ502802

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-RIPK2 (410-460) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IP, WB
Host/Source	Rabbit
Reactivity	Bovine, Human, Pig, Monkey

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.55 µg/µl
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000 IP: 1:150
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	8767
Gene Symbol	RIPK2
Uniprot ID	RIPK2_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 410-460 on human RIPK2 protein.
Immunogen Region	410-460
Specificity	
Immunogen Sequence	