

**Anti-RIPK2 antibody (410-460) {FITC} (STJ502803)**

STJ502803

**GENERAL INFORMATION**

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Rabbit polyclonal antibody anti-RIPK2 (410-460) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
<b>Applications</b>	ELISA, IP, WB
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	Bovine, Human, Pig, Monkey

**PRODUCT PROPERTIES**

<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	0.55 µg/µl
<b>Conjugation</b>	FITC
<b>Purification</b>	Affinity Purified
<b>Dilution Range</b>	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000 IP: 1:150
<b>Formulation</b>	
<b>Isotype</b>	IgG
<b>Storage Instruction</b>	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

**TARGET INFORMATION**

<b>Gene ID</b>	<a href="#">8767</a>
<b>Gene Symbol</b>	<a href="#">RIPK2</a>
<b>Uniprot ID</b>	<a href="#">RIPK2_HUMAN</a>
<b>Immunogen</b>	Synthetic peptide taken within amino acid region 410-460 on human RIPK2 protein.
<b>Immunogen Region</b>	410-460
<b>Specificity</b>	
<b>Immunogen Sequence</b>	