

Anti-RIPK3 antibody (461-511) {FITC} (STJ502807)

STJ502807

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-RIPK3 (461-511) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IP, WB
Host/Source	Rabbit
Reactivity	Bovine, Human, Pig, Monkey

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.66-0.68 µg/µl
Conjugation	FITC
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000 IP: 1:200
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	11035
Gene Symbol	RIPK3
Uniprot ID	RIPK3_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 461-511 on human Receptor-interacting Serine/Threonine-protein kinase 3 protein.
Immunogen Region	461-511
Specificity	
Immunogen Sequence	