

Anti-DSTYK antibody (200-250) {FITC} (STJ502815)

STJ502815

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-DSTYK (200-250) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IP, WB
Host/Source	Rabbit
Reactivity	Bovine, Human, Monkey, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.55 µg/µl
Conjugation	FITC
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000 IP: 1:150
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	25778
Gene Symbol	DSTYK
Uniprot ID	DUSTY_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 200-250 on human Dual serine/threonine and tyrosine protein kinase protein (RIPK5).
Immunogen Region	200-250
Specificity	
Immunogen Sequence	