

Anti-RPE65 antibody {FITC} (STJ502833)

STJ502833

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-RPE65 is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IP, WB
Host/Source	Rabbit
Reactivity	Bovine, Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.51-0.63 µg/µl
Conjugation	FITC
Purification	Affinity Purified
Dilution Range	WB: 1:500 IP: 1:200
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	19892
Gene Symbol	Rpe65
Uniprot ID	RPE65_MOUSE
Immunogen	Synthetic peptide corresponding to unique amino acid sequence on mouse RPE65 protein.
Immunogen Region	
Specificity	
Immunogen Sequence	