

Anti-RUNX1 antibody (1-80) {Biotin} (STJ502835)

STJ502835

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-RUNX1 (1-80) is suitable for use in ELISA and Western Blot research applications.
Applications	ELISA, WB
Host/Source	Rabbit
Reactivity	Bovine, Human, Monkey

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.51-0.53 µg/µl
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	WB: 1:10, 000 ELISA: 1:10, 000
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	861
Gene Symbol	RUNX1
Uniprot ID	RUNX1_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 1-80 on human RUNX1 protein.
Immunogen Region	1-80
Specificity	
Immunogen Sequence	