

Anti-RUNX3 antibody (200-250) {FITC} (STJ502845)

STJ502845

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-RUNX3 (200-250) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IP, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.52-0.53 µg/µl
Conjugation	FITC
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000 IP: 1:150
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	864
Gene Symbol	RUNX3
Uniprot ID	RUNX3_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 200-250 on human Runt-related transcription factor 3 protein.
Immunogen Region	200-250
Specificity	
Immunogen Sequence	