

## Anti-EDG8 antibody (380-400) {FITC} (STJ502857)

STJ502857

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Rabbit polyclonal antibody anti-EDG8 (380-400) is suitable for use in ELISA and Western Blot research applications.
<b>Applications</b>	ELISA, WB
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	Human, Mouse, Rat

### PRODUCT PROPERTIES

<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	0.52-0.63 µg/µl
<b>Conjugation</b>	FITC
<b>Purification</b>	Affinity Purified
<b>Dilution</b>	WB: 1:500
<b>Range</b>	DB: 1:10, 000 ELISA: 1:10, 000
<b>Formulation</b>	
<b>Isotype</b>	IgG
<b>Storage</b>	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
<b>Instruction</b>	

### TARGET INFORMATION

<b>Gene ID</b>	<a href="#">60399</a>
<b>Gene Symbol</b>	<a href="#">S1pr5</a>
<b>Uniprot ID</b>	<a href="#">S1PR5_RAT</a>
<b>Immunogen</b>	Synthetic peptide taken within amino acid region 380-400 of rat S1P5R/EDG8 protein. Does not cross react to other Spingosine 1-phosphate receptor subtypes. Labels a strong band of EDG8 of 48-50 kDa. Also labels a 54 kDa band, the identity of this is n
<b>Immunogen Region</b>	380-400
<b>Specificity</b>	
<b>Immunogen Sequence</b>	