

Anti-Adenylate Cyclase 10 antibody (200-250) {Biotin} (STJ502868)

STJ502868

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Adenylate Cyclase 10 (200-250) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IP, WB
Host/Source	Rabbit
Reactivity	Gorilla, Human, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.55-0.75 µg/µl
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000 IP: 1:150
Formulation	
Isotype	IgG
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	271639
Gene Symbol	Adcy10
Uniprot ID	ADCYA_MOUSE
Immunogen	Synthetic peptide taken within amino acid region 200-250 on mouse Adenylate Cyclase type 10 protein.
Immunogen	200-250
Region	
Specificity	
Immunogen	
Sequence	