

Anti-SEPT1 antibody (300-350) {Biotin} (STJ502915)

STJ502915

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-SEPT1 (300-350) is suitable for use in ELISA and Western Blot research applications.
Applications	ELISA, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.63-0.70 µg/µl
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1731
Gene Symbol	SEPTIN1
Uniprot ID	SEPT1_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 300-350 on SEPT1 protein.
Immunogen Region	300-350
Specificity	
Immunogen Sequence	