

Anti-SEPT2 antibody (270-320) {Biotin} (STJ502917)

STJ502917

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-SEPT2 (270-320) is suitable for use in ELISA and Western Blot research applications.
Applications	ELISA, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.64-0.72 µg/µl
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000 IP: 1:200 IHC: 1:50-1:150
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	4735
Gene Symbol	SEPTIN2
Uniprot ID	SEPT2_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 270-320 on human Septin-2 protein.
Immunogen Region	270-320
Specificity	
Immunogen Sequence	