

Anti-SEPT9 antibody (1-60) {FITC} (STJ502932)

STJ502932

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-SEPT9 (1-60) is suitable for use in ELISA and Western Blot research applications.
Applications	ELISA, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.65-0.69 µg/µl
Conjugation	FITC
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000 IP: 1:200 IHC: 1:50-1:100
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	10801
Gene Symbol	SEPTIN9
Uniprot ID	SEPT9_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 1-60 on SEPT9 protein.
Immunogen Region	1-60
Specificity	
Immunogen Sequence	