

Anti-SERINC4 antibody (350-400) {Biotin} (STJ502937)

STJ502937

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-SERINC4 (350-400) is suitable for use in ELISA, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IHC, IP, WB
Host/Source	Rabbit
Reactivity	Human, Monkey

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5 µg/µl
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	WB: 1:500 ELISA: 1:10, 000 IP: 1:200 IHC: 1:100
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	619189
Gene Symbol	SERINC4
Uniprot ID	SERC4_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 350-400 on human Serine incorporator 4 protein.
Immunogen Region	350-400
Specificity	
Immunogen Sequence	