

Anti-SESN3 antibody (1-50) {FITC} (STJ502959)

STJ502959

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-SESN3 (1-50) is suitable for use in ELISA and Western Blot research applications.
Applications	ELISA, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.68-0.80 µg/µl
Conjugation	FITC
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	143686
Gene Symbol	SESN3
Uniprot ID	SESN3_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 1-50 on SESN3 protein.
Immunogen Region	1-50
Specificity	
Immunogen Sequence	