

## Anti-Soluble Guanylyl Cyclase Alpha 2 antibody (680-730) {FITC} (STJ502974)

STJ502974

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Rabbit polyclonal antibody anti-Soluble Guanylyl Cyclase Alpha 2 (680-730) is suitable for use in ELISA, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
<b>Applications</b>	ELISA, IHC, IP, WB
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	Human, Mouse, Rat

### PRODUCT PROPERTIES

<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	0.7 µg/µl
<b>Conjugation</b>	FITC
<b>Purification</b>	Affinity Purified
<b>Dilution Range</b>	WB: 1:750 ELISA: 1:10, 000 IP: 1:200 IHC: 1:50-1:150
<b>Formulation</b>	
<b>Isotype</b>	IgG
<b>Storage Instruction</b>	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

### TARGET INFORMATION

<b>Gene ID</b>	<a href="#">66012</a>
<b>Gene Symbol</b>	<a href="#">Gucy1a2</a>
<b>Uniprot ID</b>	<a href="#">GCYA2_RAT</a>
<b>Immunogen</b>	Synthetic peptide taken within amino acid region 680-730 on rat guanylate cyclase 1, soluble, alpha 2 protein.
<b>Immunogen Region</b>	680-730
<b>Specificity</b>	
<b>Immunogen Sequence</b>	