

Anti-SGK2 antibody (377-427) {Biotin} (STJ502982)

STJ502982

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Chicken polyclonal antibody anti-SGK2 (377-427) is suitable for use in ELISA, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IHC, IP, WB
Host/Source	Chicken
Reactivity	Human, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.6 µg/µl
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000 IP: 1:200 IHC: 1:100 IF: 1:100
Formulation	
Isotype	IgY
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	10110
Gene Symbol	SGK2
Uniprot ID	SGK2_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 377-427 on human Sgk2 protein.
Immunogen Region	377-427
Specificity	
Immunogen Sequence	