

Anti-C8orf44-SGK3 antibody (440-490) (STJ502984)

GENERAL INFORMATION

Product Type Primary antibodies

Short Chicken polyclonal antibody anti-C8orf44-SGK3 (440-490) is suitable for use in Confocal Microscopy, ELISA, Immunocytochemistry,

Description Immunofluorescence, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.

Applications CM, ELISA, IC, IF, IHC, IP, WB

Host/Source Chicken

Reactivity Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality Polyclonal

Clone ID
Concentration 0.6 µg/µl Conjugation Unconjugated Purification Affinity Purified

Dilution WB: 1:500 Range DB: 1:10, 000 ELISA: 1:10, 000

IP: 1:200 IHC: 1:100 ICC: 1:100 IF: 1:100 CFM: 1:100

Formulation

Isotype IgY

Storage Store at-20°C for long term storage. Avoid freeze-thaw cycles.

Instruction

TARGET INFORMATION

Gene ID 100533105/23678

Gene Symbol SGK3

Uniprot ID SGK3_HUMAN

Immunogen Synthetic peptide taken within amino acid region 440-490 on human Sgk3 protein.

Immunogen 440-490 Region Specificity

Immunogen Sequence

