

Anti-SHANK3 antibody (1600-1650) {Biotin} (STJ502997)

STJ502997

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-SHANK3 (1600-1650) is suitable for use in ELISA, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IHC, IP, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5 µg/µl
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	WB: 1:500 ELISA: 1:10, 000 IP: 1:150 IHC: 1:50-1:100
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	58234
Gene Symbol	Shank3
Uniprot ID	SHAN3_MOUSE
Immunogen	Synthetic peptide taken within amino acid region 1600-1650 on mouse SH3 and multiple ankyrin repeat domains protein 3.
Immunogen Region	1600-1650
Specificity	
Immunogen Sequence	