

Anti-Shibire antibody (800-850) {Biotin} (STJ503009)

STJ503009

GENERAL INFORMATION

| | |
|--------------------------|--|
| Product Type | Primary antibodies |
| Short Description | Rabbit polyclonal antibody anti-Shibire (800-850) is suitable for use in ELISA and Western Blot research applications. |
| Applications | ELISA, WB |
| Host/Source | Rabbit |
| Reactivity | D.melanogaster |

PRODUCT PROPERTIES

| | |
|----------------------------|---|
| Clonality | Polyclonal |
| Clone ID | |
| Concentration | 0.50 µg/µl |
| Conjugation | Biotin |
| Purification | Affinity Purified |
| Dilution Range | WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000 IP: 1:200 |
| Formulation | |
| Isotype | IgG |
| Storage Instruction | Store at -20°C for long term storage. Avoid freeze-thaw cycles. |

TARGET INFORMATION

| | |
|---------------------------|--|
| Gene ID | |
| Gene Symbol | |
| Uniprot ID | |
| Immunogen | Synthetic peptide taken within amino acid region 800-850 on D. melanogaster Shibire protein. |
| Immunogen Region | 800-850 |
| Specificity | |
| Immunogen Sequence | |