

Anti-SIGLEC10 antibody (650-697) {Biotin} (STJ503027)

STJ503027

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-SIGLEC10 (650-697) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IP, WB
Host/Source	Rabbit
Reactivity	Human, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.56-0.60 µg/µl
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	WB: 1:500 ELISA: 1:10, 000 IP: 1:200
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	89790
Gene Symbol	SIGLEC10
Uniprot ID	SIG10_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 650-697 on human SIGLEC10 protein.
Immunogen Region	650-697
Specificity	
Immunogen Sequence	