

Anti-SIRT1 antibody (25-75) (STJ503029)

STJ503029

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-SIRT1 (25-75) is suitable for use in ELISA, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IHC, IP, WB
Host/Source	Rabbit
Reactivity	Chicken, Human, Monkey, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.69-0.75 µg/µl
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000 IP: 1:200 IHC: 1:100
Formulation	
Isotype	IgG
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	93759
Gene Symbol	Sirt1
Uniprot ID	SIR1_MOUSE
Immunogen	Synthetic peptide taken within amino acid 25-75 on rat Sirtuin 1 protein.
Immunogen Region	25-75
Specificity	
Immunogen Sequence	



Immunohistochemistry of rat tissue with Anti-SIRT1 primary antibody (STJ503029) at 1:100 dilution in Buffer. Section was treated with DAB and Haematoxylin stain.