

Anti-Solute carrier family 35 member E4 antibody (170-200) {FITC} (STJ503055)

STJ503055

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Solute carrier family 35 member E4 (170-200) is suitable for use in ELISA, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IHC, IP, WB
Host/Source	Rabbit
Reactivity	Human, Monkey, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.67-0.98 µg/µl
Conjugation	FITC
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000 IP: 1:200 IHC: 1:100
Formulation	
Isotype	IgG
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	266687
Gene Symbol	Slc35e4
Uniprot ID	S35E4_RAT
Immunogen	Synthetic peptide taken within amino acid region 170-200 on Slc35E4 protein.
Immunogen Region	170-200
Specificity	
Immunogen Sequence	