

Anti-SLC9A3 antibody (790-835) {FITC} (STJ503058)

STJ503058

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-SLC9A3 (790-835) is suitable for use in ELISA and Western Blot research applications.
Applications	ELISA, WB
Host/Source	Rabbit
Reactivity	Human, Monkey, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.68 µg/µl
Conjugation	FITC
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	6550
Gene Symbol	SLC9A3
Uniprot ID	SL9A3_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 790-835 on Slc9A3 protein.
Immunogen Region	790-835
Specificity	
Immunogen Sequence	