

Anti-SMAD4 antibody (131-181) {FITC} (STJ503064)

STJ503064

GENERAL INFORMATION

| | |
|--------------------------|---|
| Product Type | Primary antibodies |
| Short Description | Rabbit polyclonal antibody anti-SMAD4 (131-181) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications. |
| Applications | ELISA, IP, WB |
| Host/Source | Rabbit |
| Reactivity | Human, Mouse, Rat |

PRODUCT PROPERTIES

| | |
|----------------------------|--|
| Clonality | Polyclonal |
| Clone ID | |
| Concentration | 0.6 µg/µl |
| Conjugation | FITC |
| Purification | Affinity Purified |
| Dilution Range | WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000 IP: 1:200 IHC: 1:50-1:150 |
| Formulation | |
| Isotype | IgG |
| Storage Instruction | Store at -20°C for long term storage. Avoid freeze-thaw cycles. |

TARGET INFORMATION

| | |
|---------------------------|--|
| Gene ID | 4089 |
| Gene Symbol | SMAD4 |
| Uniprot ID | SMAD4_HUMAN |
| Immunogen | Synthetic peptide taken within amino acid region 131-181 on human SMAD4 protein. |
| Immunogen Region | 131-181 |
| Specificity | |
| Immunogen Sequence | |