

Anti-NAPA antibody (195-245) {FITC} (STJ503076)

STJ503076

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-NAPA (195-245) is suitable for use in ELISA, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IHC, IP, WB
Host/Source	Rabbit
Reactivity	Human, Monkey, Mouse, Pig, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.56 µg/µl
Conjugation	FITC
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000 IP: 1:150 IHC: 1:50-1:150
Formulation	
Isotype	IgG
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	8775
Gene Symbol	NAPA
Uniprot ID	SNAA_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 195-245 on human alpha-soluble NSF attachment protein.
Immunogen Region	195-245
Specificity	
Immunogen Sequence	