

## Anti-SOX10 antibody (10-60) {Biotin} (STJ503088)

STJ503088

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Rabbit polyclonal antibody anti-SOX10 (10-60) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
<b>Applications</b>	ELISA, IP, WB
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	Human, Monkey, Mouse, Rat

### PRODUCT PROPERTIES

<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	0.55-0.75 µg/µl
<b>Conjugation</b>	Biotin
<b>Purification</b>	Affinity Purified
<b>Dilution Range</b>	WB: 1:500 DB: 1:10,000 ELISA: 1:10,000 IP: 1:150 IHC: 1:50-1:150
<b>Formulation</b>	
<b>Isotype</b>	IgG
<b>Storage Instruction</b>	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

### TARGET INFORMATION

<b>Gene ID</b>	<a href="#">6663</a>
<b>Gene Symbol</b>	<a href="#">SOX10</a>
<b>Uniprot ID</b>	<a href="#">SOX10_HUMAN</a>
<b>Immunogen</b>	Synthetic peptide taken within amino acid region 10-60 on human Transcription Factor SOX-10 protein.
<b>Immunogen Region</b>	10-60
<b>Specificity</b>	
<b>Immunogen Sequence</b>	