

Anti-SOX10 antibody (400-450) {FITC} (STJ503091)

STJ503091

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-SOX10 (400-450) is suitable for use in ELISA, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IHC, IP, WB
Host/Source	Rabbit
Reactivity	Human, Monkey

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.55-0.75 µg/µl
Conjugation	FITC
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000 IP: 1:150 IHC: 1:50-1:100
Formulation	
Isotype	IgG
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	6663
Gene Symbol	SOX10
Uniprot ID	SOX10_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 400-450 on human Transcription Factor SOX-10 protein.
Immunogen Region	400-450
Specificity	
Immunogen Sequence	