

## Anti-LYZ antibody (20-70) (STJ503098)

STJ503098

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Rabbit polyclonal antibody anti-LYZ (20-70) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
<b>Applications</b>	ELISA, IP, WB
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	Human, Mouse, Rat

### PRODUCT PROPERTIES

<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	0.57 µg/µl
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	Affinity Purified
<b>Dilution Range</b>	WB: 1:500 ELISA: 1:10, 000 IP: 1:200
<b>Formulation</b>	
<b>Isotype</b>	IgG
<b>Storage Instruction</b>	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

### TARGET INFORMATION

<b>Gene ID</b>	4069
<b>Gene Symbol</b>	LYZ
<b>Uniprot ID</b>	LYSC_HUMAN
<b>Immunogen</b>	Synthetic peptides taken within amino acid region 20-70 on human sperm acrosome membrane-associated protein 3 isoform 2.
<b>Immunogen Region</b>	20-70
<b>Specificity</b>	This antibody will label all isoforms of sperm acrosome membrane-associated protein 3.
<b>Immunogen Sequence</b>	