

Anti-SPACA3 antibody (50-100) {FITC} (STJ503106)

STJ503106

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-SPACA3 (50-100) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IP, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.57 µg/µl
Conjugation	FITC
Purification	Affinity Purified
Dilution Range	WB: 1:500 ELISA: 1:10, 000 IP: 1:200 IHC: 1:50
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	124912
Gene Symbol	SPACA3
Uniprot ID	SACA3_HUMAN
Immunogen	Synthetic peptides taken within amino acid region 50-100 and 80-130 on human sperm acrosome membrane-associated protein 3 isoform 2.
Immunogen Region	50-100
Specificity	This antibody will label all isoforms of sperm acrosome membrane-associated protein 3.
Immunogen Sequence	