

Anti-SPACA3 antibody (50-100) (STJ503107)

STJ503107

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-SPACA3 (50-100) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IP, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.57 µg/µl
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB: 1:500 ELISA: 1:10, 000 IP: 1:200 IHC: 1:50
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	124912
Gene Symbol	SPACA3
Uniprot ID	SACA3_HUMAN
Immunogen	Synthetic peptides taken within amino acid region 50-100 and 80-130 on human sperm acrosome membrane-associated protein 3 isoform 2.
Immunogen Region	50-100
Specificity	This antibody will label all isoforms of sperm acrosome membrane-associated protein 3.
Immunogen Sequence	



Immunohistochemistry of rat tissue with Anti-SPACA3 primary antibody (STJ503107) at 1:50 dilution in Buffer. Section was treated with DAB (brown) staining and Haematoxylin QS (blue) counterstain. Magnification is at 5X.