

Anti-SREBF1 antibody (950-1000) {Biotin} (STJ503130)

STJ503130

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-SREBF1 (950-1000) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IP, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.65 µg/µl
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000 IP: 1:150
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	6720
Gene Symbol	SREBF1
Uniprot ID	SRBP1_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 950-1000 on human Sterol regulatory element binding transcription factor 1 protein.
Immunogen Region	950-1000
Specificity	
Immunogen Sequence	