

Anti-SREBF2 antibody (640-690) {FITC} (STJ503133)

STJ503133

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-SREBF2 (640-690) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IP, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.64-0.72 µg/µl
Conjugation	FITC
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000 IP: 1:150
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	6721
Gene Symbol	SREBF2
Uniprot ID	SRBP2_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 640-690 on human Sterol regulatory element binding protein-2.
Immunogen Region	640-690
Specificity	
Immunogen Sequence	