

Anti-STRADA antibody (381-431) (STJ503156)

STJ503156

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-STRADA (381-431) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IP, WB
Host/Source	Rabbit
Reactivity	Human, Monkey

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.64-0.72 µg/µl
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000 IP: 1:150
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	92335
Gene Symbol	STRADA
Uniprot ID	STRAA_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 381-431 on human STE-20 related kinase adapter protein alpha.
Immunogen Region	381-431
Specificity	
Immunogen Sequence	