

## Anti-STRADA antibody (381-431) {FITC} (STJ503160)

STJ503160

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Rabbit polyclonal antibody anti-STRADA (381-431) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
<b>Applications</b>	ELISA, IP, WB
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	Human, Monkey

### PRODUCT PROPERTIES

<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	0.64-0.72 µg/µl
<b>Conjugation</b>	FITC
<b>Purification</b>	Affinity Purified
<b>Dilution Range</b>	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000 IP: 1:150
<b>Formulation</b>	
<b>Isotype</b>	IgG
<b>Storage Instruction</b>	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

### TARGET INFORMATION

<b>Gene ID</b>	92335
<b>Gene Symbol</b>	STRADA
<b>Uniprot ID</b>	STRAA_HUMAN
<b>Immunogen</b>	Synthetic peptide taken within amino acid region 381-431 on human STE-20 related kinase adapter protein alpha.
<b>Immunogen Region</b>	381-431
<b>Specificity</b>	
<b>Immunogen Sequence</b>	