

Anti-BIRC5 antibody (STJ503164)

STJ503164

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-BIRC5 is suitable for use in ELISA and Western Blot research applications.
Applications	ELISA, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.50-1.50 µg/µl
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	332
Gene Symbol	BIRC5
Uniprot ID	BIRC5_HUMAN
Immunogen	Synthetic peptide corresponding to unique amino acid sequences on human Survivin protein.
Immunogen Region	
Specificity	
Immunogen Sequence	