

Anti-Syncytin 2 antibody (380-430) (STJ503173)

STJ503173

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Syncytin 2 (380-430) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IP, WB
Host/Source	Rabbit
Reactivity	Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5 µg/µl
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB: 1:500 ELISA: 1:50, 000 IP: 1:200
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	239167
Gene Symbol	Synb
Uniprot ID	SYNB_MOUSE
Immunogen	Synthetic peptide taken within amino acid region 380-430 on human Envelope glycoprotein syncytin-B protein.
Immunogen Region	380-430
Specificity	
Immunogen Sequence	



Western blot with Anti-Syncytin 2 Antibody (STJ503173) on rat: 1) Cortex, 2) Midbrain, 3) Cerebellum, 4) Hippocampus, 5) Kidney, 6) Adrenal gland, 7) Thymus, 8) Trachea, 9) Ovary, 10) Bone marrow, 11) Testis. Antibody at 1:250 dilution in buffer.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081