

Anti-ERVFRD-1 antibody (80-130) {Biotin} (STJ503174)

STJ503174

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-ERVFRD-1 (80-130) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IP, WB
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5 µg/µl
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	WB: 1:500 ELISA: 1:50, 000 IP: 1:200
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	405754
Gene Symbol	ERVFRD-1
Uniprot ID	SYCY2_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 80-130 on human Syncytin-2 protein.
Immunogen Region	80-130
Specificity	
Immunogen Sequence	