

Anti-SYK antibody (100-150) {Biotin} (STJ503183)

STJ503183

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-SYK (100-150) is suitable for use in Confocal Microscopy, ELISA, Immunocytochemistry, Immunofluorescence, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
Applications	CM, ELISA, ICC, IF, IHC, IP, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.66 µg/µl
Conjugation	Biotin
Purification	Affinity Purified
Dilution	WB: 1:500
Range	ELISA: 1:10, 000 IP: 1:200 IHC: 1:50-1:200 ICC: 1:50-1:200 IF: 1:50-1:200
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	6850
Gene Symbol	SYK
Uniprot ID	KSYK_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 100-150 on human Syk protein.
Immunogen Region	100-150
Specificity	
Immunogen Sequence	