

Anti-SYP antibody (1-100) (STJ503185)

GENERAL INFORMATION

Product Type Primary antibodies

Short Rabbit polyclonal antibody anti-SYP (1-100) is suitable for use in Confocal Microscopy, ELISA, Immunocytochemistry,

Description Immunofluorescence, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.

Applications CM, ELISA, ICC, IF, IHC, IP, WB

Host/Source Rabbit

Reactivity Chicken, Human, Monkey, Mouse, Rat

PRODUCT PROPERTIES

Clonality Polyclonal

Clone ID Concentration 0.60-0.62 µg/µl Conjugation Unconjugated Purification Affinity Purified **Dilution** WB: 1:1, 000 **Range** DB: 1:10, 000 ELISA: 1:10, 000

> IP: 1:125 IHC: 1:100 ICC: 1:100 IF: 1:100 CFM: 1:100

Formulation

Isotype IgG

Storage Store at-20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID 6855

Gene Symbol SYP

Uniprot ID SYPH_HUMAN

Immunogen Synthetic peptide taken within amino acid region 1-100 on human Synaptophysin protein.

Immunogen 1-100

Region

Specificity This antibody does not with cross-react other synaptic vesicle proteins such as synaptotagmin, GAP45, and SNAP25 proteins. Immunogen

Sequence

