

Anti-SDC1 antibody (250-350) {Biotin} (STJ503195)

STJ503195

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-SDC1 (250-350) is suitable for use in ELISA and Western Blot research applications.
Applications	ELISA, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.55-0.61 µg/µl
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	6382
Gene Symbol	SDC1
Uniprot ID	SDC1_HUMAN
Immunogen	C-epitope synthetic peptide taken within amino acid region 250-350 on human Syndecan-1 protein.
Immunogen Region	250-350
Specificity	
Immunogen Sequence	