

Anti-SYT1 antibody (1-60) {Biotin} (STJ503198)

STJ503198

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-SYT1 (1-60) is suitable for use in ELISA and Western Blot research applications.
Applications	ELISA, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.64-0.72 µg/µl
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	6857
Gene Symbol	SYT1
Uniprot ID	SYT1_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 1-60 on SYT1 protein.
Immunogen Region	1-60
Specificity	
Immunogen Sequence	