

Anti-TDRD1 antibody (570-605) {Biotin} (STJ503225)

STJ503225

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-TDRD1 (570-605) is suitable for use in ELISA, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IHC, IP, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.52-0.53 µg/µl
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000
Formulation	
Isotype	IgG
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	56165
Gene Symbol	TDRD1
Uniprot ID	TDRD1_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 570-605 on TDRD1 protein.
Immunogen Region	570-605
Specificity	
Immunogen Sequence	