

Anti-TDRKH antibody (480-515) {Biotin} (STJ503228)

STJ503228

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-TDRKH (480-515) is suitable for use in ELISA and Western Blot research applications.
Applications	ELISA, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.56-0.58 µg/µl
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000 IP: 1:150 IHC: 1:50-1:150
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	11022
Gene Symbol	TDRKH
Uniprot ID	TDRKH_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 480-515 on TDRD2 protein.
Immunogen Region	480-515
Specificity	
Immunogen Sequence	