

Anti-TGM1 antibody {FITC} (STJ503242)

GENERAL INFORMATION

Product Type Primary antibodies

Short Rabbit polyclonal antibody anti-TGM1 is suitable for use in Confocal Microscopy, ELISA, Immunocytochemistry,

Description Immunohistochemistry, Immunofluorescence and Western Blot research applications.

Applications CM, ELISA, ICC, IHC, IF, WB

Host/Source Rabbit

Reactivity Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality Polyclonal

Clone ID Concentration 0.61-0.65 µg/µl

Conjugation FITC

Purification Affinity Purified **Dilution** WB: 1:450 **Range** DB: 1:30, 000

ELISA: 1:30, 000 IHC: 1:250 ICC: 1:250 IF: 1:250 CFM: 1:250

Formulation

Isotype IgG

Storage Store at-20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID 74145

Gene Symbol F13a1 Uniprot ID F13A_MOUSE

Immunogen Synthetic peptide corresponding to unique amino acid sequence on Coagulation factor XIIIA chain precursor

Immunogen Region Specificity Immunogen Sequence