

Anti-TGM1 antibody {FITC} (STJ503244)

STJ503244

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-TGM1 is suitable for use in Confocal Microscopy, ELISA, Immunocytochemistry, Immunohistochemistry, Immunofluorescence and Western Blot research applications.
Applications	CM, ELISA, ICC, IHC, IF, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.61-0.65 µg/µl
Conjugation	FITC
Purification	Affinity Purified
Dilution Range	WB: 1:450 DB: 1:30, 000 ELISA: 1:30, 000 IHC: 1:250 ICC: 1:250 IF: 1:250 CFM: 1:250
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	74145
Gene Symbol	F13a1
Uniprot ID	F13A_MOUSE
Immunogen	Synthetic peptide corresponding to unique amino acid sequence on Coagulation factor XIIIa chain precursor
Immunogen Region	
Immunogen Specificity	
Immunogen Sequence	