

Anti-TEK antibody (540-590) (STJ503258)

STJ503258

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Chicken polyclonal antibody anti-TEK (540-590) is suitable for use in Dot Blot, ELISA, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
Applications	DB, ELISA, IHC, IP, WB
Host/Source	Chicken
Reactivity	Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.68-0.70 µg/µl
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:4, 000 ELISA: 1:4, 000 IP: 1:150 IHC: 1:50-1:200
Formulation	
Isotype	IgY
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	7010
Gene Symbol	TEK
Uniprot ID	TIE2_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 540-590 on human Tyrosine-protein kinase receptor Tie-1 isoform 2 protein.
Immunogen Region	540-590
Specificity	
Immunogen Sequence	



Western blot with Anti-TEK Antibody (STJ503258) with peptide on 1) Mouse brain, 2) TIE2 peptide. Antibody at 1:4000 dilution in buffer.