

Anti-Thickveins antibody {FITC} (STJ503271)

STJ503271

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Thickveins is suitable for use in ELISA and Western Blot research applications.
Applications	ELISA, WB
Host/Source	Rabbit
Reactivity	D.melanogaster

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.70 µg/µl
Conjugation	FITC
Purification	Affinity Purified
Dilution Range	DB-1:10, 000 ELISA-1:10, 000 IP-1:200 WB-1:500
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Immunogenic peptide was covalently modified to achieve the desired antigenicity.
Immunogen Region	
Specificity	
Immunogen Sequence	