

Anti-TLR1 antibody (370-420) (STJ503276)

STJ503276

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-TLR1 (370-420) is suitable for use in ELISA, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IHC, IP, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.55-0.67 µg/µl
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB: 1:500 ELISA: 1:20, 000 IP: 1:200 IHC: 1:100 ICC: 1:100 IF: 1:100
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	7096
Gene Symbol	TLR1
Uniprot ID	TLR1_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 370-420 on human Toll-like receptor 1.
Immunogen Region	370-420
Specificity	
Immunogen Sequence	



Immunohistochemistry of human thymus with Anti-TLR1 primary antibody (STJ503276) at 1:700 dilution in Buffer.