

## Anti-TLR5 antibody (195-225) (STJ503279)

STJ503279

### GENERAL INFORMATION

|                          |  |
|--------------------------|--|
| <b>Product Type</b>      | Primary antibodies   |
| <b>Short Description</b> | Rabbit polyclonal antibody anti-TLR5 (195-225) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications. |
| <b>Applications</b>      | ELISA, IP, WB  |
| <b>Host/Source</b>       | Rabbit   |
| <b>Reactivity</b>        | Human  |

### PRODUCT PROPERTIES

|                            |   |
|----------------------------|---|
| <b>Clonality</b>           | Polyclonal  |
| <b>Clone ID</b>            |   |
| <b>Concentration</b>       | 0.75 µg/µl  |
| <b>Conjugation</b>         | Unconjugated  |
| <b>Purification</b>        | Affinity Purified   |
| <b>Dilution Range</b>      | WB: 1:500<br>DB: 1:20, 000<br>ELISA: 1:20, 000<br>IP: 1:200     |
| <b>Formulation</b>         |   |
| <b>Isotype</b>             | IgG   |
| <b>Storage Instruction</b> | Store at -20°C for long term storage. Avoid freeze-thaw cycles. |

### TARGET INFORMATION

|                           |   |
|---------------------------|---|
| <b>Gene ID</b>            | 7100  |
| <b>Gene Symbol</b>        | TLR5  |
| <b>Uniprot ID</b>         | TLR5_HUMAN  |
| <b>Immunogen</b>          | Synthetic peptide taken within amino acid region 195-225 on human TLR5 protein. |
| <b>Immunogen Region</b>   | 195-225   |
| <b>Specificity</b>        |   |
| <b>Immunogen Sequence</b> |   |